

Configuration des 2 BD LSRFortessa X-20 :

Laser	P M T	Miroir dichroïque	Filtre BP	Option de fluorochromes
Bleu 488 nm	A	655 LP	695/30	PerCP-Cy5.5, 7-AAD, PerCP
	B	505 LP	530/30	FITC, GFP, CFSE, Alexa488, Live/Dead Green, YFP
Violet 405 nm	A	750LP	780/60	BV785
	B	685 LP	710/50	BV711, efluor 700
	C	635 LP	670/30	BV650
	D	600LP	610/20	eFluor 605, BV605
	E	505LP	525/50	AmCyan, Alexa Fluor 430, V500, BV510, V500, Pacific Orange, LiveDead Aqua
	F	420LP	450/50	Pacific Blue, Alexa Fluor 405, V450, BV421, Cascade Blue, eFluor450, DAPI, CFP...
Red 633 nm	A	705LP	780/60	APC-CY7, APC-H7, APC-Alexa750, Alexa750, APC-eFluor780
	B	690LP	730/45	Alexa-Fluor700, APC-Cy5.5, APC-Alexa700, eFluor710, AlexaFluor680
	C	665LP	670/30	APC, Alexa647, APC, Alexa633
Yellow/Green 561nm	A	750LP	780/60	PE-CY7
	B	685LP	710/50	PE-Alexa 700, PE-Cy5-5
	C	600LP	610/20	PE-Alexa 594, PE-Alexa 610, PE-CF594, PE-mCherry, PE-Texas Red, PI
	D	EMPTY	586/15	DsRed, PE, Alexa 546
UV 355 nm	A	690LP	740/35	BUV737
	B	450LP	515/30	BUV496
	C	370LP	379/28	BUV395